

折页：
尺寸：492*165mm
157铜版纸

<div data-bbox="856 397 1068 439" data-label="Image"> </div> <div data-bbox="856 449 1068 519" data-label="Image"> </div> <div data-bbox="856 526 1068 582" data-label="Text"> <p>Cincom www.cincomhealth.com Email: service@cincomhealth.com</p> </div> <div data-bbox="789 624 1145 689" data-label="Text"> <p>SHENZHEN CINCOM E-COMMERCE CO., LTD. Address: B806, building 11, gongyuandaxiaocu, Longxiangdadao ongchengjiedao, longgangqu, shenzhen, guangdong, China E-mail: service@cincomhealth.com Website: www.cincomhealth.com</p> </div> <div data-bbox="789 689 1024 726" data-label="Image"> </div>	<div data-bbox="1246 282 1384 301" data-label="Section-Header"> <p>TWO MODES OPTIONAL</p> </div> <div data-bbox="1246 301 1602 375" data-label="Text"> <p>1. Mode A: Synchronous Mode When you press the power button, the machine will automatically work in mode A. It will provide kneading and shiatsu massage in the low air pressure intensity, you can change the air pressure intensity (L1-L2-L3) and heat (H1-H2-OF) according to your own preferences.</p> </div> <div data-bbox="1246 375 1525 412" data-label="Text"> <p>Feature: In this mode, the air pressure and kneading function will work simultaneously in cycle.</p> </div> <div data-bbox="1246 412 1602 467" data-label="Text"> <p>2. Mode B: Composite Mode This mode provides a combination massage and will also run in the low air pressure intensity. You can change the air pressure intensity (L1-L2-L3) and heat (H1-H2-OF) according to your own preferences.</p> </div> <div data-bbox="1246 467 1602 522" data-label="Text"> <p>Feature: This mode will provide a different massage method. Massage functions will work in turns: air pressure massage mode-kneading massage mode-combination massage (like mode A), running in cycle.</p> </div> <div data-bbox="1246 522 1602 559" data-label="Text"> <p>Note: The foot shiatsu massager will be automatically shut off in 10/20/30 minutes (3 massage timings optional), you can restart it anytime if you need.</p> </div> <div data-bbox="1246 587 1491 606" data-label="Section-Header"> <p>PRODUCT MAINTENANCE AND ATTENTIONS</p> </div> <div data-bbox="1246 615 1602 670" data-label="List-Group"> <ul style="list-style-type: none"> • Please disconnect the cable from the power socket after use. • Ensure the device is disconnected before cleaning. • Do not clean the device with a damp cloth, brushes or hard materials which might damage the same. • Do not use the device for more than 40 minutes a day. </div> <div data-bbox="1246 726 1344 744" data-label="Section-Header"> <p>SPECIFICATION</p> </div> <div data-bbox="1246 744 1454 800" data-label="List-Group"> <ul style="list-style-type: none"> • Product Name: Shiatsu Foot Massager • Model No: CM-002F • Voltage: DC12V/3.0A • Power: 40W </div>	<div data-bbox="1720 282 1841 301" data-label="Section-Header"> <p>TROUBLESHOOTING</p> </div> <div data-bbox="1720 301 2076 569" data-label="Table"> <table border="1"> <thead> <tr> <th>PROBLEM</th> <th>POSSIBLE</th> <th>CAUSE/SOLUTION</th> </tr> </thead> <tbody> <tr> <td rowspan="3">Does not work</td> <td>The power supply is not on</td> <td>Turn on the power supply</td> </tr> <tr> <td>Power plug is poorly connected</td> <td>Check the connection</td> </tr> <tr> <td>The power plug has been damaged</td> <td>Inspect the plug for any damage. Consult the professional service team for repair.</td> </tr> <tr> <td rowspan="2">A sudden interruption of the power</td> <td>Power button has not been pressed</td> <td>Press the on/off button</td> </tr> <tr> <td>Power off</td> <td>Check whether the power supply interface is connected well</td> </tr> <tr> <td></td> <td>Message timing is over</td> <td>Press the on/off button to restart if you need more massage</td> </tr> </tbody> </table> </div> <div data-bbox="1720 578 1797 596" data-label="Section-Header"> <p>CONTACT US</p> </div> <div data-bbox="1720 596 2076 624" data-label="Text"> <p>We provide 2-year warranty and long-term after-sale service. If you have any problem during using of this product, please feel free to contact us.</p> </div> <div data-bbox="1720 624 1898 643" data-label="Text"> <p>E-mail: service@cincomhealth.com</p> </div> <div data-bbox="1720 643 2076 670" data-label="Text"> <p>Note: Please write the order number together with the problems you meet in the mail</p> </div>	PROBLEM	POSSIBLE	CAUSE/SOLUTION	Does not work	The power supply is not on	Turn on the power supply	Power plug is poorly connected	Check the connection	The power plug has been damaged	Inspect the plug for any damage. Consult the professional service team for repair.	A sudden interruption of the power	Power button has not been pressed	Press the on/off button	Power off	Check whether the power supply interface is connected well		Message timing is over	Press the on/off button to restart if you need more massage	<div data-bbox="2412 282 2553 319" data-label="Image"> </div> <div data-bbox="2284 328 2553 365" data-label="Section-Header"> <p>Operating Instruction</p> </div> <div data-bbox="2284 375 2553 430" data-label="Section-Header"> <p>SHIATSU FOOT MASSAGER CM-002F</p> </div> <div data-bbox="2268 449 2479 726" data-label="Image"> </div> <div data-bbox="2217 781 2520 837" data-label="Text"> <p>Thank you for choosing the CINCOM Foot massager. Please read this manual carefully and learn how to handle the device in the correct manner before using it. Pay special attention to the safety instructions and keep this manual for future reference.</p> </div>
PROBLEM	POSSIBLE	CAUSE/SOLUTION																			
Does not work	The power supply is not on	Turn on the power supply																			
	Power plug is poorly connected	Check the connection																			
	The power plug has been damaged	Inspect the plug for any damage. Consult the professional service team for repair.																			
A sudden interruption of the power	Power button has not been pressed	Press the on/off button																			
	Power off	Check whether the power supply interface is connected well																			
	Message timing is over	Press the on/off button to restart if you need more massage																			
<div data-bbox="799 892 873 911" data-label="Section-Header"> <p>CONTENTS</p> </div> <div data-bbox="799 920 1152 1049" data-label="Table-Of-Contents"> <p>Introduction — 2 Safety precautions — 2 Structure of the device — 4 Operation instructions — 5 Product maintenance and attentions — 6 Specification — 6 Troubleshooting — 7 Contact us — 7</p> </div> <div data-bbox="799 1058 1152 1114" data-label="Text"> <p>Thank you for choosing the CINCOM Foot massager. Please read this manual carefully and learn how to handle the device in the correct manner before using it. Pay special attention to the safety instructions and keep this manual for future reference.</p> </div> <div data-bbox="799 1123 890 1142" data-label="Section-Header"> <p>INTRODUCTION</p> </div> <div data-bbox="799 1151 1152 1317" data-label="List-Group"> <ul style="list-style-type: none"> • Nice appearance, portable and convenient. • Combines with shiatsu massage, kneading massage, air pressure massage and heating function together. • Kneading & shiatsu massage from toes to heels to cover all acupuncture points. • Air pressure massage makes you more relaxed. • 2 Massage modes: user mode and auto mode available • 3 Levels of air pressure intensities, three options for massage timing, two working mode available. • 2 Levels of heating temperatures adjustment: low, high temperature. • 10/20/30 minutes auto shut-off function. • Powered by a DC12V adapter, safe and reliable. </div> <div data-bbox="799 1363 934 1382" data-label="Section-Header"> <p>SAFETY PRECAUTIONS</p> </div> <div data-bbox="799 1391 1152 1419" data-label="Text"> <p>The product should be used in correct manners, as explained in this manual, to prevent physical and material damage.</p> </div>	<div data-bbox="1246 892 1327 911" data-label="Section-Header"> <p>WARNINGS</p> </div> <div data-bbox="1246 911 1602 948" data-label="Text"> <p>People with any of the following conditions or those who are receiving medical treatments should consult the doctor prior to using the shiatsu foot massager:</p> </div> <div data-bbox="1246 948 1602 1142" data-label="List-Group"> <ul style="list-style-type: none"> • Being pregnant or at puerperium • Having a fever • The doctor has recommended you rest or you are not feeling well. • Convalescing on the doctors orders. • Suffering from a chronic illness or cancer. • Having cerebral lesions (clots, lacerations, traumatism, etc.) • Recovering from an operation to the upper part of your body, buttocks, thighs or calves. • Suffering from a cardiac disease, heart-related problems or use an electronic medical device. • Using a pacemaker. • Injured or have a skin disease. • Having broken or malformed bones. • This device is not suitable for children, the disabled or invalids. • The use of the product is not recommended for people with psychic problems. </div> <div data-bbox="1246 1151 1293 1169" data-label="Section-Header"> <p>NOTE:</p> </div> <div data-bbox="1246 1169 1602 1206" data-label="Text"> <p>If you feel the Air Pressure intensity too big to make you uncomfortable during use, we advise you to press the "Deflate" button for 5 seconds to relieve the discomfort.</p> </div> <div data-bbox="1246 1215 1327 1234" data-label="Section-Header"> <p>WHILE IN USE</p> </div> <div data-bbox="1246 1234 1602 1391" data-label="List-Group"> <ul style="list-style-type: none"> • Do not put the device in a place which are damp or extremely dusty, as this may result in mechanical failure. • Do not place the device in sources of heat or in the sunlight. • Place the reflexotherapy device on a flat surface. • Do not use the device together with other therapeutic equipment or electrical medical apparatus and similar. • Please avoid any liquid splashing on the equipment when using the equipment • Do not use the device on any wet part of the body, as this may result in mechanical failure or even an electric shock. • Your body should be relaxed during treatment. • Do not use the device immediately after meals, wait for at least 1 hour. • Daily use should be limited to a maximum of 30 minutes. </div> <div data-bbox="1246 1391 1327 1410" data-label="Section-Header"> <p>Attention:</p> </div> <div data-bbox="1246 1410 1602 1428" data-label="Text"> <p>Check the voltage is suitable and never remove the cable directly.</p> </div>	<div data-bbox="1720 892 1898 911" data-label="Section-Header"> <p>WHEN THE DEVICE IS NOT IN USE</p> </div> <div data-bbox="1720 911 2076 948" data-label="List-Group"> <ul style="list-style-type: none"> • Ensure the device is turned off after use. • Unplug the device from the mains after use and before cleaning. • Store the device in a clean and dry location. </div> <div data-bbox="1720 966 1898 984" data-label="Section-Header"> <p>THE DEVICE SHOULD NOT BE USED</p> </div> <div data-bbox="1720 984 2076 1058" data-label="List-Group"> <ul style="list-style-type: none"> • If it falls into water; the device should not be placed in or come into contact with water or other liquids. • If you feel strange or extremely uncomfortable when using the device, turn it off immediately and call your doctor. • If there is electrical storm. • If the plug or cable is damaged or malfunctioning. </div> <div data-bbox="1720 1086 1898 1105" data-label="Section-Header"> <p>STRUCTURE OF THE DEVICE</p> </div> <div data-bbox="1720 1132 2076 1410" data-label="Diagram"> </div>	<div data-bbox="2177 892 2335 911" data-label="Section-Header"> <p>OPERATION INSTRUCTION</p> </div> <div data-bbox="2177 911 2530 957" data-label="Text"> <p>Use the power adapter to connect the foot massager and the power socket, then it will make a "di" sound, and all indicators on the panel flash. The foot massager is ready for working now.</p> </div> <div data-bbox="2177 984 2352 1003" data-label="Section-Header"> <p>FUNCTIONS OF CONTROL PANEL</p> </div> <div data-bbox="2177 1012 2496 1142" data-label="Diagram"> </div> <div data-bbox="2177 1206 2530 1410" data-label="List-Group"> <ol style="list-style-type: none"> 1. Press the Power button to start massage (the indicator shows fixed light), the machine will run in the lowest intensity. Mode A and 20-minute massage timing by default. Note: The heating function is closed by default, you can press the heat button to start it during use if you need. 2. Press the Time button to select massage timing: 10, 20 or 30 minutes. 3. Press the Air Pressure button to select air pressure intensity, "L1-L2-L3" (low-medium-high) can be seen through the display screen on the panel. 4. Press the Heat button to select low or high temperature or cancel heating; "H1-H2-OF" will be showed on the display screen on the panel. 5. Press the Mode B button to end the massage, the indicator will flash. </div>																		